

POU domain, class 2, transcription factor 2

Catalog No: #AP86707



Package Size: #AP86707-1 50ug #AP86707-2 100ug #AP86707-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	POU domain, class 2, transcription factor 2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Lymphoid-restricted immunoglobulin octamer-binding protein NF-A2, Octamer-binding protein 2, Oct-2, Octamer-binding transcription factor 2, OTF-2, Oct2, Otf2
Accession No.	Uniprot ID:Q00196
Uniprot	Q00196
GeneID	18987;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: Pou2f2

Note: This product is for in vitro research use only