

Nuclear receptor subfamily 4 group A member 1

Catalog No: #AP86708



Package Size: #AP86708-1 50ug #AP86708-2 100ug #AP86708-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Nuclear receptor subfamily 4 group A member 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Early response protein NAK1,Nuclear hormone receptor NUR/77,Nur77,Orphan nuclear receptor HMR,Orphan nuclear receptor TR3,ST-59,Testicular receptor 3,GFRP1,HMR,NAK1
Accession No.	Uniprot ID:P22736
Uniprot	P22736
GeneID	3164;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:NR4A1

Note: This product is for in vitro research use only