

Neurotrophin-4

Catalog No: #AP86709



Package Size: #AP86709-1 50ug #AP86709-2 100ug #AP86709-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

| | |
|--------------------|--|
| Product Name | Neurotrophin-4 |
| Brief Description | Recombinant Protein |
| Host Species | E.coli |
| Purification | Greater than 90% by SDS-PAGE |
| Species Reactivity | Hu |
| Other Names | NT-4,Neurotrophin-5,NT-5,Neutrophic factor 4,NTF5 |
| Accession No. | Uniprot ID:P34130 |
| Uniprot | P34130 |
| GeneID | 4909; |
| Tag Info | His |
| Storage | Store at -20°C. (Avoid repeated freezing and thawing.) |

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:NTF4

Note: This product is for in vitro research use only