

## Serum paraoxonase/arylesterase 2

Catalog No: #AP86714



Package Size: #AP86714-1 50ug #AP86714-2 100ug #AP86714-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Serum paraoxonase/arylesterase 2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Bovine
Other Names	PON 2,Aromatic esterase 2,A-esterase 2,Serum aryldialkylphosphatase 2
Accession No.	Uniprot ID:Q58DS7
Uniprot	Q58DS7
GeneID	281417;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:PON2

Note: This product is for in vitro research use only