

Serum paraoxonase/arylesterase 2

Catalog No: #AP86715



Package Size: #AP86715-1 50ug #AP86715-2 100ug #AP86715-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Serum paraoxonase/arylesterase 2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	PON 2,Aromatic esterase 2,A-esterase 2,Serum aryldialkylphosphatase 2
Accession No.	Uniprot ID:Q6AXM8
Uniprot	Q6AXM8
GeneID	296851;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Pon2

Note: This product is for in vitro research use only