

Prostaglandin-H2 D-isomerase

Catalog No: #AP86719



Package Size: #AP86719-1 50ug #AP86719-2 100ug #AP86719-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Prostaglandin-H2 D-isomerase
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Dog
Other Names	Glutathione-independent PGD synthase,Lipocalin-type prostaglandin-D synthase,Prostaglandin-D2 synthase,PGD2 synthase,PGDS,PGDS2
Accession No.	Uniprot ID:Q9XS65
Uniprot	Q9XS65
GeneID	403740;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:PTGDS

Note: This product is for in vitro research use only