

Persephin

Catalog No: #AP86721



Package Size: #AP86721-1 50ug #AP86721-2 100ug #AP86721-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Persephin
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	PSP
Accession No.	Uniprot ID:O60542
Uniprot	O60542
GeneID	5623;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: PSPN

Note: This product is for in vitro research use only