

1-phosphatidylinositol 4,5-bisphosphate phosphodiesterase beta-1

Catalog No: #AP86725

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Package Size: #AP86725-1 50ug #AP86725-2 100ug #AP86725-3 1mg

Description

Product Name	1-phosphatidylinositol 4,5-bisphosphate phosphodiesterase beta-1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Bovine
Other Names	PLC-154,Phosphoinositide phospholipase C-beta-1,Phospholipase C-beta-1,PLC-beta-1
Accession No.	Uniprot ID:P10894
Uniprot	P10894
GeneID	287026;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:PLCB1

Note: This product is for in vitro research use only