

[Pyruvate dehydrogenase (acetyl-transferring)] kinase isozyme 1, mitochondrial

Catalog No: #AP86726

Orders: order@signalwayantibody.com

Package Size: #AP86726-1 50ug #AP86726-2 100ug #AP86726-3 1mg

Support: tech@signalwayantibody.com

Description

Product Name	[Pyruvate dehydrogenase (acetyl-transferring)] kinase isozyme 1, mitochondrial
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Pyruvate dehydrogenase kinase isoform 1,PDH kinase 1
Accession No.	Uniprot ID:Q8BFP9
Uniprot	Q8BFP9
GeneID	228026;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Pdk1

Note: This product is for in vitro research use only