

## Xaa-Pro dipeptidase

Catalog No: #AP86729



Package Size: #AP86729-1 50ug #AP86729-2 100ug #AP86729-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Xaa-Pro dipeptidase
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	X-Pro dipeptidase,Imidodipeptidase,Peptidase 4,Peptidase D,Proline dipeptidase,Prolidase,Pep4
Accession No.	Uniprot ID:Q11136
Uniprot	Q11136
GeneID	18624;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:Pepd

Note: This product is for in vitro research use only