

Pituitary adenylate cyclase-activating polypeptide

Catalog No: #AP86732



Package Size: #AP86732-1 50ug #AP86732-2 100ug #AP86732-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Pituitary adenylate cyclase-activating polypeptide
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Pig
Other Names	PACAP
Accession No.	Uniprot ID:P41535
Uniprot	P41535
GeneID	414283;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: ADCYAP1

Note: This product is for in vitro research use only