

Phosphatidylinositol-glycan-specific phospholipase D

Catalog No: #AP86737



Package Size: #AP86737-1 50ug #AP86737-2 100ug #AP86737-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Phosphatidylinositol-glycan-specific phospholipase D
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	PI-G PLD, Glycoprotein phospholipase D, Glycosyl-phosphatidylinositol-specific phospholipase D, GPI-PLD, PIGPLD1
Accession No.	Uniprot ID:P80108
Uniprot	P80108
GeneID	2822;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: GPLD1

Note: This product is for in vitro research use only