

Peroxiredoxin-6

Catalog No: #AP86738



Package Size: #AP86738-1 50ug #AP86738-2 100ug #AP86738-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Peroxiredoxin-6
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	1-Cys peroxiredoxin,1-Cys PRX,Acidic calcium-independent phospholipase A2,aiPLA2,Antioxidant protein 2,Non-selenium glutathione peroxidase,NSGPx,Aop2,Ltw4,Prdx5
Accession No.	Uniprot ID:O08709
Uniprot	O08709
GeneID	11758;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Prdx6

Note: This product is for in vitro research use only