

## Glycodelin

Catalog No: #AP86740



Package Size: #AP86740-1 50ug #AP86740-2 100ug #AP86740-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Glycodelin
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	GD,Placental protein 14,PP14,Pregnancy-associated endometrial alpha-2 globulin,PAEG,Progestagen-associated endometrial protein,Progesterone-associated endometrial protein
Accession No.	Uniprot ID:P09466
Uniprot	P09466
GeneID	5047;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:PAEP

Note: This product is for in vitro research use only