

Serine/threonine-protein kinase PAK 1

Catalog No: #AP86741



Package Size: #AP86741-1 50ug #AP86741-2 100ug #AP86741-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Serine/threonine-protein kinase PAK 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Alpha-PAK,CDC42/RAC effector kinase PAK-A,p21-activated kinase 1,PAK-1,p65-PAK,Paka
Accession No.	Uniprot ID:O88643
Uniprot	O88643
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Pak1

Note: This product is for in vitro research use only