

## Prostaglandin D2 receptor

Catalog No: #AP86742



Package Size: #AP86742-1 50ug #AP86742-2 100ug #AP86742-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Prostaglandin D2 receptor
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	PGD receptor, Prostanoid DP receptor
Accession No.	Uniprot ID:P70263
Uniprot	P70263
GeneID	19214;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer: 50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene\_full\_name: Ptgdr

Note: This product is for in vitro research use only