

Protein disulfide-isomerase

Catalog No: #AP86743



Package Size: #AP86743-1 50ug #AP86743-2 100ug #AP86743-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Protein disulfide-isomerase
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	PDI, Cellular thyroid hormone-binding protein, Prolyl 4-hydroxylase subunit beta,p55,ERBA2L,PDI,PDIA1,PO4DB
Accession No.	Uniprot ID:P07237
Uniprot	P07237
GeneID	5034;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: P4HB

Note: This product is for in vitro research use only