

## Tartrate-resistant acid phosphatase type 5

Catalog No: #AP86746



Package Size: #AP86746-1 50ug #AP86746-2 100ug #AP86746-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Tartrate-resistant acid phosphatase type 5
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	TR-AP, Tartrate-resistant acid ATPase, TrATPase, Type 5 acid phosphatase, T5ap, Trap
Accession No.	Uniprot ID:Q05117
Uniprot	Q05117
GeneID	11433;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer: 50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene\_full\_name: Acp5

Note: This product is for in vitro research use only