

Alkaline phosphatase, placental type

Catalog No: #AP86747



Package Size: #AP86747-1 50ug #AP86747-2 100ug #AP86747-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Alkaline phosphatase, placental type
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Alkaline phosphatase Regan isozyme,Placental alkaline phosphatase 1,PLAP-1,PLAP
Accession No.	Uniprot ID:P05187
Uniprot	P05187
GeneID	250;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:ALPP

Note: This product is for in vitro research use only