

Prolactin-inducible protein

Catalog No: #AP86748



Package Size: #AP86748-1 50ug #AP86748-2 100ug #AP86748-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Prolactin-inducible protein
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Gross cystic disease fluid protein 15,GCDFP-15,Prolactin-induced protein,Secretory actin-binding protein,SABP,gp17,GCDFP15,GPIP4
Accession No.	Uniprot ID:P12273
Uniprot	P12273
GeneID	5304;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:PIP

Note: This product is for in vitro research use only