

Plasminogen activator inhibitor 2, macrophage

Catalog No: #AP86750



Package Size: #AP86750-1 50ug #AP86750-2 100ug #AP86750-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Plasminogen activator inhibitor 2, macrophage
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	PAI-2,Serpin B2,Pai2,Planh2
Accession No.	Uniprot ID:P12388
Uniprot	P12388
GeneID	18788;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Serpinb2

Note: This product is for in vitro research use only