

Proteinase-activated receptor 1

Catalog No: #AP86751



Package Size: #AP86751-1 50ug #AP86751-2 100ug #AP86751-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Proteinase-activated receptor 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Bovine
Other Names	PAR-1,Thrombin receptor,PAR1
Accession No.	Uniprot ID:A7YY44
Uniprot	A7YY44
GeneID	526585;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:F2R

Note: This product is for in vitro research use only