

1-phosphatidylinositol 4,5-bisphosphate phosphodiesterase gamma-2

Catalog No: #AP86758

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Package Size: #AP86758-1 50ug #AP86758-2 100ug #AP86758-3 1mg

Description

Product Name	1-phosphatidylinositol 4,5-bisphosphate phosphodiesterase gamma-2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Phosphoinositide phospholipase C-gamma-2,Phospholipase C-gamma-2,PLC-gamma-2
Accession No.	Uniprot ID:Q8CIH5
Uniprot	Q8CIH5
GeneID	234779;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Plcg2

Note: This product is for in vitro research use only