

Pituitary adenylate cyclase-activating polypeptide

Catalog No: #AP86764



Package Size: #AP86764-1 50ug #AP86764-2 100ug #AP86764-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Pituitary adenylate cyclase-activating polypeptide
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	PACAP,Pacap
Accession No.	Uniprot ID:O70176
Uniprot	O70176
GeneID	11516;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Adcyap1

Note: This product is for in vitro research use only