

Phospholipase A2, membrane associated

Catalog No: #AP86767



Package Size: #AP86767-1 50ug #AP86767-2 100ug #AP86767-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Phospholipase A2, membrane associated
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Bovine
Other Names	GIIC sPLA2, Group IIA phospholipase A2, Phosphatidylcholine 2-acylhydrolase 2A
Accession No.	Uniprot ID:Q56JZ2
Uniprot	Q56JZ2
GeneID	507201;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: PLA2G2A

Note: This product is for in vitro research use only