

Parathyroid hormone/parathyroid hormone-related peptide receptor

Catalog No: #AP86770

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Package Size: #AP86770-1 50ug #AP86770-2 100ug #AP86770-3 1mg

Description

Product Name	Parathyroid hormone/parathyroid hormone-related peptide receptor
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	PTH/PTHrP type I receptor,PTH/PTHr receptor,Parathyroid hormone 1 receptor,PTH1 receptor,Pthr,Pthr1
Accession No.	Uniprot ID:P25961
Uniprot	P25961
GeneID	56813;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Pth1r

Note: This product is for in vitro research use only