

Glycogen phosphorylase, brain form

Catalog No: #AP86773



Package Size: #AP86773-1 50ug #AP86773-2 100ug #AP86773-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Glycogen phosphorylase, brain form
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Accession No.	Uniprot ID:Q8CI94
Uniprot	Q8CI94
GeneID	110078;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Pygb

Note: This product is for in vitro research use only