

Glutaminyl-peptide cyclotransferase-like protein

Catalog No: #AP86774



Package Size: #AP86774-1 50ug #AP86774-2 100ug #AP86774-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Glutaminyl-peptide cyclotransferase-like protein
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Golgi-resident glutaminyl-peptide cyclotransferase,isoQC,gQC
Accession No.	Uniprot ID:Q9NXS2
Uniprot	Q9NXS2
GeneID	54814;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:QPCTL

Note: This product is for in vitro research use only