

Oligodendrocyte transcription factor 2

Catalog No: #AP86776



Package Size: #AP86776-1 50ug #AP86776-2 100ug #AP86776-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Oligodendrocyte transcription factor 2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Oligo2,Class B basic helix-loop-helix protein 1,bHLHb1,Class E basic helix-loop-helix protein 19,bHLHe19,Protein kinase C-binding protein 2,Protein kinase C-binding protein RACK17,BHLHB1,BHLHE19,PRKCBP2,RACK17
Accession No.	Uniprot ID:Q13516
Uniprot	Q13516
GeneID	10215;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:OLIG2

Note: This product is for in vitro research use only