

Osteomodulin

Catalog No: #AP86777



Package Size: #AP86777-1 50ug #AP86777-2 100ug #AP86777-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Osteomodulin
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	Keratan sulfate proteoglycan osteomodulin,KSPG osteomodulin,Osteoadherin,OSAD
Accession No.	Uniprot ID:Q9Z1S7
Uniprot	Q9Z1S7
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH2PO4, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Omd

Note: This product is for in vitro research use only