

Pleckstrin homology-like domain family A member 2

Catalog No: #AP86780



Package Size: #AP86780-1 50ug #AP86780-2 100ug #AP86780-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Pleckstrin homology-like domain family A member 2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Beckwith-Wiedemann syndrome chromosomal region 1 candidate gene C protein,Imprinted in placenta and liver protein,Tumor-suppressing STF cDNA 3 protein,Tumor-suppressing subchromosomal transferable fragment candidate gene 3 protein,p17-Beckwith-Wiedemann region 1 C,p17-BWR1C,BWR1C,HLDA2,IPL,TSSC3
Accession No.	Uniprot ID:Q53GA4
Uniprot	Q53GA4
GeneID	7262;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:PHLDA2

Note: This product is for in vitro research use only