

## Pulmonary surfactant-associated protein C

Catalog No: #AP86782



Package Size: #AP86782-1 50ug #AP86782-2 100ug #AP86782-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Pulmonary surfactant-associated protein C
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	SP-C,Pulmonary surfactant-associated proteolipid SPL(Val),Sftp2
Accession No.	Uniprot ID:P11685
Uniprot	P11685
GeneID	50683;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:Sftpc

Note: This product is for in vitro research use only