

1-phosphatidylinositol 4,5-bisphosphate phosphodiesterase gamma-1

Catalog No: #AP86785

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Package Size: #AP86785-1 50ug #AP86785-2 100ug #AP86785-3 1mg

Description

Product Name	1-phosphatidylinositol 4,5-bisphosphate phosphodiesterase gamma-1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	Phosphoinositide phospholipase C-gamma-1, Phospholipase C-gamma-1, PLC-gamma-1
Accession No.	Uniprot ID:P10686
Uniprot	P10686
GeneID	25738;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: Plcg1

Note: This product is for in vitro research use only