

## N-fatty-acyl-amino acid synthase/hydrolase PM20D1

Catalog No: #AP86787



Package Size: #AP86787-1 50ug #AP86787-2 100ug #AP86787-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	N-fatty-acyl-amino acid synthase/hydrolase PM20D1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Chicken
Other Names	Peptidase M20 domain-containing protein 1,RCJMB04_8d17
Accession No.	Uniprot ID:Q5ZL18
Uniprot	Q5ZL18
GeneID	419838;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer: 50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene\_full\_name: PM20D1

Note: This product is for in vitro research use only