

Prolyl 4-hydroxylase subunit alpha-3

Catalog No: #AP86788



Package Size: #AP86788-1 50ug #AP86788-2 100ug #AP86788-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Prolyl 4-hydroxylase subunit alpha-3
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Bovine
Other Names	4-PH alpha-3,Procollagen-proline,2-oxoglutarate-4-dioxygenase subunit alpha-3
Accession No.	Uniprot ID:Q75UG4
Uniprot	Q75UG4
GeneID	414348;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:P4HA3

Note: This product is for in vitro research use only