

# Patatin-like phospholipase domain-containing protein 2

Catalog No: #AP86789

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Package Size: #AP86789-1 50ug #AP86789-2 100ug #AP86789-3 1mg

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Patatin-like phospholipase domain-containing protein 2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Adipose triglyceride lipase, Calcium-independent phospholipase A2, Desnutrin, IPLA2-zeta, Pigment epithelium-derived factor, TTS2.2, Transport-secretion protein 2, TTS2, ATGL, FP17548
Accession No.	Uniprot ID:Q96AD5
Uniprot	Q96AD5
GeneID	57104;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer: 50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene\_full\_name: PNPLA2

Note: This product is for in vitro research use only