

## N-acetylmuramoyl-L-alanine amidase

Catalog No: #AP86790



Package Size: #AP86790-1 50ug #AP86790-2 100ug #AP86790-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	N-acetylmuramoyl-L-alanine amidase
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Peptidoglycan recognition protein 2,Peptidoglycan recognition protein long,PGRP-L,TagL,Pglyrpl,Pgrpl
Accession No.	Uniprot ID:Q8VCS0
Uniprot	Q8VCS0
GeneID	57757;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:Pglyrp2

Note: This product is for in vitro research use only