

# Protein-L-isoaspartate(D-aspartate) O-methyltransferase

Catalog No: #AP86791

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

Package Size: #AP86791-1 50ug #AP86791-2 100ug #AP86791-3 1mg

## Description

Product Name	Protein-L-isoaspartate(D-aspartate) O-methyltransferase
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	PIMT,L-isoaspartyl protein carboxyl methyltransferase,Protein L-isoaspartyl/D-aspartyl methyltransferase,Protein-beta-aspartate methyltransferase
Accession No.	Uniprot ID:P23506
Uniprot	P23506
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:Pcmt1

Note: This product is for in vitro research use only