

Orexin-A

Catalog No: #AP86792



Package Size: #AP86792-1 50ug #AP86792-2 100ug #AP86792-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Orexin-A
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Bovine
Other Names	Hypocretin-1,Hcrt1,OX,PPOX
Accession No.	Uniprot ID:P56717
Uniprot	P56717
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH2PO4, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:HCRT

Note: This product is for in vitro research use only