

N-fatty-acyl-amino acid synthase/hydrolase PM20D1

Catalog No: #AP86794



Package Size: #AP86794-1 50ug #AP86794-2 100ug #AP86794-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	N-fatty-acyl-amino acid synthase/hydrolase PM20D1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Peptidase M20 domain-containing protein 1
Accession No.	Uniprot ID:Q8C165
Uniprot	Q8C165
GeneID	212933;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: Pm20d1

Note: This product is for in vitro research use only