

cAMP-specific 3',5'-cyclic phosphodiesterase 4B

Catalog No: #AP86798



Package Size: #AP86798-1 50ug #AP86798-2 100ug #AP86798-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	cAMP-specific 3',5'-cyclic phosphodiesterase 4B
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	DPDE4,PDE32,DPDE4
Accession No.	Uniprot ID:Q07343
Uniprot	Q07343
GeneID	5142;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:PDE4B

Note: This product is for in vitro research use only