

Protein disulfide-isomerase A4

Catalog No: #AP86804



Package Size: #AP86804-1 50ug #AP86804-2 100ug #AP86804-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Protein disulfide-isomerase A4
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	Calcium-binding protein 2,CaBP2,Endoplasmic reticulum resident protein 70,ER protein 70,Endoplasmic reticulum resident protein 72,ER protein 72,Cabp2,Erp70
Accession No.	Uniprot ID:P38659
Uniprot	P38659
GeneID	116598;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Pdia4

Note: This product is for in vitro research use only