

Platelet factor 4

Catalog No: #AP86805



Package Size: #AP86805-1 50ug #AP86805-2 100ug #AP86805-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Platelet factor 4
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	PF-4,C-X-C motif chemokine 4,Cxcl4,Scyb4
Accession No.	Uniprot ID:Q9Z126
Uniprot	Q9Z126
GeneID	56744;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name: Pf4

Note: This product is for in vitro research use only