

Patatin-like phospholipase domain-containing protein

2

Catalog No: #AP86807

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Package Size: #AP86807-1 50ug #AP86807-2 100ug #AP86807-3 1mg

Description

Product Name	Patatin-like phospholipase domain-containing protein 2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Bovine
Other Names	Adipose triglyceride lipase,ATGL
Accession No.	Uniprot ID:Q2K118
Uniprot	Q2K118
GeneID	508493;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:PNPLA2

Note: This product is for in vitro research use only