

Serum paraoxonase/arylesterase 1

Catalog No: #AP86809



Package Size: #AP86809-1 50ug #AP86809-2 100ug #AP86809-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Serum paraoxonase/arylesterase 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	PON 1,Aromatic esterase 1,A-esterase 1,Serum aryldialkylphosphatase 1,Pon
Accession No.	Uniprot ID:P55159
Uniprot	P55159
GeneID	84024;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Pon1

Note: This product is for in vitro research use only