

Protein phosphatase 1 regulatory subunit 14B

Catalog No: #AP86810



Package Size: #AP86810-1 50ug #AP86810-2 100ug #AP86810-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Protein phosphatase 1 regulatory subunit 14B
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Phosphatase holoenzyme inhibitor 1,PHI-1,Phospholipase C-beta-3 neighbouring gene protein,Png
Accession No.	Uniprot ID:Q62084
Uniprot	Q62084
GeneID	18938;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Ppp1r14b

Note: This product is for in vitro research use only