

Neurofilament heavy polypeptide

Catalog No: #AP86812



Package Size: #AP86812-1 50ug #AP86812-2 100ug #AP86812-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Neurofilament heavy polypeptide
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	NF-H,200 kDa neurofilament protein,Neurofilament triplet H protein,KIAA0845,NFH
Accession No.	Uniprot ID:P12036
Uniprot	P12036
GeneID	4744;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:NEFH

Note: This product is for in vitro research use only