

Neurofilament medium polypeptide

Catalog No: #AP86814



Package Size: #AP86814-1 50ug #AP86814-2 100ug #AP86814-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

| | |
|--------------------|--|
| Product Name | Neurofilament medium polypeptide |
| Brief Description | Recombinant Protein |
| Host Species | E.coli |
| Purification | Greater than 90% by SDS-PAGE |
| Species Reactivity | Bovine |
| Other Names | NF-M, 160 kDa neurofilament protein, Neurofilament 3, Neurofilament triplet M protein, NEF3, NFM |
| Accession No. | Uniprot ID:O77788 |
| Uniprot | O77788 |
| Tag Info | His |
| Storage | Store at -20°C. (Avoid repeated freezing and thawing.) |

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: NEFM

Note: This product is for in vitro research use only