

## Oxidized low-density lipoprotein receptor 1

Catalog No: #AP86815



Package Size: #AP86815-1 50ug #AP86815-2 100ug #AP86815-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Oxidized low-density lipoprotein receptor 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Bovine
Other Names	Ox-LDL receptor 1, Lectin-like oxidized LDL receptor 1, LOX-1, Lectin-type oxidized LDL receptor 1, LOX1
Accession No.	Uniprot ID:P79391
Uniprot	P79391
GeneID	281368;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer: 50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene\_full\_name: OLR1

Note: This product is for in vitro research use only