

Peroxiredoxin-2

Catalog No: #AP86817



Package Size: #AP86817-1 50ug #AP86817-2 100ug #AP86817-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Peroxiredoxin-2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Natural killer cell-enhancing factor B,NKEF-B,PRP,Thiol-specific antioxidant protein,TSA,Thioredoxin peroxidase 1,Thioredoxin-dependent peroxide reductase 1,NKEFB,TDPX1
Accession No.	Uniprot ID:P32119
Uniprot	P32119
GeneID	7001;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:PRDX2

Note: This product is for in vitro research use only